



# forest dental

affordable excellence + designer friendly

## Forest 3901 Patient Chair



**Made in  
the USA**

**Designed for  
Australian  
requirements**

With over 35 years' experience building simple pneumatically controlled dental delivery systems in the USA, Forest is renowned for its quality at a competitive price.



### FEATURES

- Ultraleather asepsis upholstery standard
- Hydraulic cantilever operation
- Seat height 43cm (17") to 74cm (29")
- Emergency safety shutoff
- 30 degree rotating seat with lock
- Multifunction touchpad with 4 programmable positions, left/right mountable
- One hand adjustable double articulating headrest
- Standard twin backrest with lumbar support
- Left/right patient access
- Swing-down locking armrests

**SEE OVER FOR OPTIONS & EXTRAS**

PO Box 126 Blackheath, Beautiful Blue Mountains, NSW 2785 Australia

Phone: +612 4787 5382 • Fax: +612 4787 7469

Email: [info@forestdental.com.au](mailto:info@forestdental.com.au) • Web: [www.forestdental.com.au](http://www.forestdental.com.au)

# Forest 3901 Patient Chair

## OPTIONS

### **NARROW BACK**



Provides unrestricted patient access, perfect for pediatric or orthodontic dentistry.

### **NARROW BACK WITH SLINGS**



Slings provide an alternative for additional patient elbow support and create a barrier between patient and practitioner.

### **PEDO BOOSTER SEAT**



Accommodates smaller patients, ideal for children.

### **ADJUSTABLE IV BOARD**

Double articulating scissor arm accommodates any size patient.

### **SOFT-TOUCH BACK MASSAGE**

Available for standard backrest

Standard back  
28" wide



Narrow back  
24" wide



Narrow back 3907-999-NB



Narrow back with slings 907-999-NS

## UPHOLSTERY OPTIONS

Forest chairs come standard with Ultraleather asepsis upholstery, in a choice of more than a dozen different colours.

Options include:

- **Plush upholstery**
- **Memory Comfort plush upholstery**  
Designed for long procedures to give greater patient comfort, Memory Comfort consists of 3 layers of foam topped with 1.5 inches of memory foam. A Memory Comfort chair features side-bolsters on the seat cushion support legs and a memory foam headrest.
- **Single or dual colours**
- **Raffia pattern fabric**
- **Custom colours**

A range of matching Doctor and Assistant stools is also available.



Pedo booster seat 3907-180



Adjustable IV board 3900-200

*To find out more about our chairs & packages, call us on (02) 4787 1455*

PO Box 126 Blackheath, Beautiful Blue Mountains, NSW 2785 Australia

Phone: +612 4787 5382 • Fax: +612 4787 7469

Email: [info@forestdental.com.au](mailto:info@forestdental.com.au) • Web: [www.forestdental.com.au](http://www.forestdental.com.au)